



ENHANCED PROTECTION WITH GREATER VISIBILITY

BIO-HYALURTMCS

Sodium Hyaluronate 2% &
Sodium Chondroitin Sulfate 2% Ophthalmic Solution



www.biotechvisioncare.com

BIO-HYALURTMCS

Sodium Hyaluronate 2% &
Sodium Chondroitin Sulfate 2% Ophthalmic Solution

Features

- **Combination of Sodium Hyaluronate + Sodium Chondroitin Sulfate**
- **Excellent protection to Corneal Endothelium**
- **Excellent Chamber maintenance & stability**
- **Very quiet eyes post-operatively**

Specifications

Raw Material	Sodium Hyaluronate content (2.0% W/V) & Sodium Chondroitin Sulfate content (2.0% W/V)
Viscosity	@ 1 Shear Rate - 40,000 ± 15,000 mPas
pH	6.8 to 7.6
Osmolarity	340 - 400 mOsm/kg
Shelf life	2 years from manufacture
Storage	Store between 2°C to 25°C. Do not Freeze. Protect from light
Application / Intended Use	Anterior segment procedures



Bio-Tech Vision Care Pvt. Ltd.

Manufacturing Facility

555 / 556 / 557, nr. Subham Tex-O-Pack,
Khatraj, dist. Gandhinagar, Gujarat, India

India Business

210, Sarthik - II, Opp. Rajpath Club,
Sarkhej-Gandhinagar highway,
Ahmedabad-380054, Gujarat, India
Ph. / Fax: 079-40038374, 32420133
E-mail: sales@biotechvisioncare.com

International Business

401, Sarthik - II, Opp. Rajpath Club,
Sarkhej-Gandhinagar highway,
Ahmedabad-380054, Gujarat, India
Ph.: +91-79-26870896/97/98
Fax: +91-79-40033767
E-mail: intbusiness@biotechvisioncare.com

www.biotechvisioncare.com